SUPERFORM TRAY SYSTEM - Reusable impression holders in polycarbonate for edentulous and partially edentulous mouths.
System specifically designed for taking high precision impressions in the production of mobile dentures.
Applying the traditional systems of elaboration of the dentures, a perfect result can be achieved using Superform anatomical impression holders, equipped with tubercles behind the upper and lower molars to facilitate correct pressure on the gums, which is very important for a perfect fit of the dentures.
Only for professional dental use.

TECHNICAL FEATURES:

Material: Reusable polycarbonate that can be sterilized in a steam autoclave at 140°C.
IMPORTANT: The result of the process of sterilization is the duty and responsibility of the professional, who should control and verify this process in conformity with the following harmonized European Standard:
EN 554 (UNI EN ISO 554) Sterilization of medical devices. Method for validation and systematic control of steam sterilization.
The use of adhesives is unnecessary thanks to the presence of three systems of mechanical retention on the surface of the impression holder:

- Dotted surface
- Stabilizing projections
- Perforations

The careful design takes into consideration all the anatomical details of the gum arch.
Each product is hand finished and individually checked.

ADVANTAGES

The anatomical shape of the dental arch combined with the provision of tubercles behind the molars ensures a technically perfect impression for the 13 different jaw sizes, which as regards jaw impression holders, these reproduce the irregularities, sliding delicately around the spaces at the rear of the gums without irritating the palate.
The handle is specially designed by permit external manipulation by the professional without interfering with the lips of the patient.
With Superform, it is possible to achieve the accuracy of individual/anatomical impression holders without the cost of the individual impression holder.
Polycarbonate is a material that makes the products light and easy to manage and at the same time highly resistant to breakage and deformation.
Impression holders of the Superform line can be used with all types of materials for making impressions, but for best results we recommend the following materials:
Chromatic alginate, Kromalgin alginate, Protesil alginate and the Protesil/Prestige line of silicon.

TECHNOLOGY

A calliper is provided for making correct measurement of the space between the corners of the mouth of the patient or between the corners and any existing teeth. This makes it possible to select the most suitable impression holder for the patient more accurately and rapidly.

PACKAGING

Superform Tray System (Cod. 101024)

Kit containing:
13 lower moulds for edentulous mouths
7 upper moulds for edentulous mouths
4 moulds for partially edentulous mouths
1 steel calliper for making measurements

Superform Tray System COLOR (Cod. 101025)

Kit containing:
13 lower moulds for edentulous mouths (orange)
7 upper moulds for edentulous mouths (green)
4 moulds for partially edentulous mouths (red)
1 steel calliper for making measurements

The Superform Tray System kits contain 24 pieces that have been statistically selected to cover the widest possible range of cases of dental configuration of patients.
In the case of Superform Color, the differentiation by colours facilitates the choice of the mould needed.

CONSERVATION

It is advisable to store the product in a dry place.
The lot number is indicated on all packages. In communicating with the supplier, it is important to mention this information so as to make it possible to identify the lot involved promptly and correctly.

PACKING

Box containing 1 Superform Kit
Carton containing 5 kits

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